

# Tourniquet

WWW.WOJRAT.PL

## Range of application


For the vascular hemorrhage on wounded limbs, the battle line emergency first aid, also used in pre-hospital emergency treatment

## Specification

Unit:mm

Specification	basic sizes			
	width	Tolerances	length	Tolerances
38×950mm	38	±2	950	±10

## Operating procedures



- A Rub Buckle
- B Rod
- C Fixed Buckle
- D Adhesive Tape
- E Fixed Slot

### Arm Bleeding One-handed Operation

- ① Take it out of the box and open the Adhesive tape, sheathe it to the Wound closed to the heart
- ② Tighten the adhesive tape, stick tight Reversely
- ③ Stuck tight along the upper arm, don't Cover the screw rods
- ④ Rotating the rod, until stopped bleeding
- ⑤ Fix the rod in slot, then the bleeding Have been controlled
- ⑥ As to the upper limb or calf, continue plus tightly and placed on the rod
- ⑦ Lock the screw rods and adhesive tape Fixed clasp

1



2



3



4



5



6




7





## Double-hand operation legs bleeding

- ① When arounding the legs, let the free end through the inside buttonhole
- ② Back cross the tape through outer buttonhole, fix the band with friction
- ③ Tighten the adhesive tape, reversely tighted
- ④ Rotating the rod, until the bleeding stopped
- ⑤ Stuck the rod with the fixed slot, the then bleeding have been controlled
- ⑥ Lock the screw rods with the fixed clasp



## Usage of friction buckle

When manually operated a tourniquet or when the tourniquet is serious polluted, use special design of friction belt buckle to fix the band.

- a. Let the adhesive tape free end through buckles inside
- b. Back through the outer buttonhole, so the adhesive tape can be firmly positioning
- c. Tighten the adhesive tape, reversely tighted

## Cautions

- a. Disposables, single use
- b. This tourniquet be used by medical first-aid trained personnel.
- c. Read the instructions before using.
- d. Medical waste disposal requirements: according to local regulations.
- e. According to medical wastes treatment, not discarded.

## Storage

- a. Under the storage conditions, the period of validity is 5 years.
- b. Generally should be stored in a temperature -5-30 degrees Celsius, relative humidity no more than 80%, non-corrosive gases and ventilated interior.